
The FREELINER® cup includes a 2.6mm peripheral rim meant for:
 - preventing the stem neck from impinging on the inner edge of the ceramic liner
 - securing the fixation of the PEXEL® and PEXEL® -E liners while preserving the integrity of the conical seating.

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The instrumentation of the FREELINER® cup allows the operator to pre-impact the liner on the table BEFORE introducing it into the patient's acetabulum. This possibility helps to control the axis of introduction of the liner into the cup, to make sure the fixation taper is clean and dry, and to easily calibrate the impaction force. This surgical protocol which is specific to the FREELINER® cup has been developed in accordance with the Safety Reminder of the ceramic material manufacturer [1] which recommends to:

- make sure that the acetabular cup and the ceramic insert are clean dry and not damaged,
- carefully assemble the components,
- confirm proper assembly and then impact,

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E v o l u t i s
 C R E A T E U R F A B R I C A N T

Acetabular cup for Ceramic or Highly-Crosslinked PEXEL® -E Polyethylene liner

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Evolutis
MOTION INSIDE

[1] Safety Reminder / BioloX inserts. Ceramtec GmbH, Medical Products Division, D-73207 Plochingen, Germany, August 2014.
http://www.ceramtec.com/files/mt_safety_reminder_inserts_en.pdf

The **FREELINER**® acetabular cup for Total Hip Arthroplasty is designed for cementless fixation to the bone, and enable 3 choices of friction bearing:

- moderately crosslinked PEXEL® polyethylene against a metal or a ceramic head,
- highly-crosslinked polyethylene with vitamin-E (PEXEL®-E) against a metal or a ceramic head,
- ceramic composite against a ceramic composite head.

The **FREELINER**® cup has been designed to ensure maximum safety during implantation and security for the patient when used with either polyethylene or ceramic liners. The cup is uncompromisingly adapted to the different characteristics of these 2 bearings:

- the inner surface of the taper is smooth and free of asperities to facilitate and secure the insertion of the ceramic insert,
- the stability of the polyethylene insert is secured by 10 anti-rotational spurs and a locking ring situated on the equator.

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Freeliner®

- Choice of 3 friction bearings
- Primary stability enhanced by the high roughness of the coating
- Osseo-integration on highly porous substrate (40 to 70%)
- Safety and reproducibility of the surgical technique

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Moderately crosslinked PEXEL® polyethylene

- Moderately crosslinked at 40kGy max.
- VacUpac vacuum packaging
- Anti-dislocation 6mm Posterior wall
- Anterior chamfer of 32°
- Inner Ø28mm

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Highly-crosslinked PEXEL®-E polyethylene with vitamin-E

- Sequentially crosslinked UHMWPE at 100kGy max.
- VacUpac vacuum packaging
- Oxydation inhibited by content of 0.1%wt of α -tocopherol (vitamin-E) blended in the material before compression molding
- with circumferential chamfer of 20°
- Inner Ø32 and 36mm

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Ceramic (Composite of Alumina and Zirconia)

- Composite of alumina and zirconia ceramic
- New SLWZero design with constant radius
- Center of rotation aligned with the equator of the liner
- Liner protected from "chipping" by the peripheral border on the cup
- Inner Ø32, 36, and 40mm

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Hemispheric cup
Primary stability through a plasma spray high rough grade porous coating
Dimensional press-fit from 1.26 to 1.66mm according to the diameter of the cup
Optional complementary fixation with screws (Ø6.0 cancellous self-tapping):

- 4 holes in the superior part of the cup on a constant bow of 120°
- angular clearance of 20° for the screws

Cementless osseo-integration coating:

- 80µm of calcium hydroxyapatite
- 400 to 600µm of porous titanium

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