


Evolutis

CREATEUR FABRICANT



		Length (mm)		Diameter (mm)				
A	Metaphysis	65		20				
		75						
		85						
		95						
B	Straight Stem	115	10	12	14	16	18	20
		145	10	12	14	16	18	20
		175	10	12	14	16	18	20
C	Curved Stem	205	10	12	14	16	18	20
D	Lateral offset	43	(with +0mm ball head)					
E	Locking Screw	25	4.5					
		30						
		35						
		40						

PRIUS Revision Femoral Implant Dimension chart